
EUROCAVE **

Articulated Presentation Kit

Ref.: AOPRESAR/POPRESAR

Your EuroCave cabinet was designed to evolve.


This articulated presentation kit, assemble onto a «Main du Sommelier» sliding shelf, enables you to combine storage, easy access to bottles and sloping presentation.




I) Place the storage shelf supports on the interior walls at the back of the cabinet.




2) Position the shelf on the supports.



3) Place the shelf guides on the inner wall of the cabinet above the sliding shelf (leave 3 spaces free between the shelf and the guide).



4) Position the cradle on the **unloaded** sliding shelf and screw on.





To check installation, pull the shelf forwards, the cradle should be guided in the sliders.

Your kit + shelf installation can hold 20 bottles in storage, 6 bottles on the sliding shelf and 6 bottles presented.









Warning: Bottles should never be in contact with the back of the cabinet.

Cabinet model	Possible Layout	Cabinet capacity
V 283 / S 283 / 5283V / 5283S	4 Kits + Step	142 bottles
E 283	2 Kits + 3 sliders + 1 storage + Compartment	155 bottles
D 283	3 Kits + slider + storage + Compartment	152 bottles
V 183 / S 183	2 Kits + 3 sliders + Step	114 bottles
E 183	I Kit + 3 sliders + I storage + Compartment	123 bottles
D 183	I Kit + 5 sliders + I storage + Compartment	136 bottles
V 083 / S 083 / 5083V / 5083S	Kit + slider + Step	58 bottles
V 292	3 Kits + 3 sliders + Step	138 bottles
5290	3 Kits + 1 slider + Step	I I 4 bottles